

# Knopfgewichte



MESSING\* feingedreht

## EINZELGEWICHTE

Nennwert	Tol.	ohne Etui	Reinraum Schraubdose	Reinraum Kunststoff-Etui	Holz-Etui Buche	Service Koffer	DAkKS Kalibrierschein
	+/- mg	Art.-Nr.	Art.-Nr.	Art.-Nr.	Art.-Nr.	Art.-Nr.	4.0000-Art.-Nr.
1 g	3	4.KMFO-110	4.KMFD-110	4.KMFP-110	4.KMFB-110		...111
2 g	4	4.KMFO-120	4.KMFD-120	4.KMFP-120	4.KMFB-120		...121
5 g	5	4.KMFO-130	4.KMFD-130	4.KMFP-130	4.KMFB-130		...131
10 g	6	4.KMFO-140	4.KMFD-140	4.KMFP-140	4.KMFB-140		...141
20 g	8	4.KMFO-150	4.KMFD-150	4.KMFP-150	4.KMFB-150		...151
50 g	10	4.KMFO-160	4.KMFD-160	4.KMFP-160	4.KMFB-160		...161
100 g	16	4.KMFO-170	4.KMFD-170		4.KMFB-170		...171
200 g	30	4.KMFO-180	4.KMFD-180		4.KMFB-180		...181
500 g	80	4.KMFO-190	4.KMFD-190		4.KMFB-190		...191
1 kg	160	4.KMFO-210	4.KMFD-210		4.KMFB-210		...211
2 kg	300	4.KMFO-220	4.KMFD-220		4.KMFB-220		...221
5 kg	800	4.KMFO-230	4.KMFD-230		4.KMFB-230		...231
10 kg	1600	4.KMFO-240	4.KMFD-240		4.KMFB-240	4.KMFS-240	...241
20 kg	3000	4.KMFO-250			4.KMFB-250	4.KMFS-250	...251

## SÄTZE

Nennwert	Anzahl	Inhalt	Reinraum Kunststoff Koffer	Holz-Etui Buche	Service Koffer	DAkKS Kalibrierschein
			Art.-Nr.	Art.-Nr.	Art.-Nr.	Art.-Nr.
1 g – 50 g	8	110 g	4.XAFP-560	4.XAFB-560	4.XAFS-560	4.0000-561
1 g – 100 g	9	210 g	4.XAFP-570	4.XAFB-570	4.XAFS-570	4.0000-571
1 g – 200 g	11	610 g	4.XAFP-580	4.XAFB-580	4.XAFS-580	4.0000-581
1 g – 500 g	12	1.110 g	4.XAFP-590	4.XAFB-590	4.XAFS-590	4.0000-591
1 g – 1 kg	13	2.110 g	4.XAFP-610	4.XAFB-610	4.XAFS-610	4.0000-611
1 g – 2 kg	15	6.110 g	4.XAFP-620	4.XAFB-620	4.XAFS-620	4.0000-621
1 g – 5 kg	16	11.110 g	4.XAFP-630	4.XAFB-630	4.XAFS-630	4.0000-631
1 g – 10 kg	17	21.110 g	4.XAFP-640	4.XAFB-640	4.XAFS-640	4.0000-641

## SÄTZE mit Milligrammstücken (Plättchenform)

Nennwert	Anzahl	Inhalt	Reinraum Kunststoff Koffer	Holz-Etui Buche	Service Koffer	DAkKS Kalibrierschein
			Art.-Nr.	Art.-Nr.	Art.-Nr.	Art.-Nr.
100 mg – 50 g	20	110,9 g	4.XMFP-760	4.XMFB-760	4.XMFS-760	4.0000-761
100 mg – 100 g	21	210,9 g	4.XMFP-770	4.XMFB-770	4.XMFS-770	4.0000-771
100 mg – 200 g	23	610,9 g	4.XMFP-780	4.XMFB-780	4.XMFS-780	4.0000-781
100 mg – 500 g	24	1.110,9 g	4.XMFP-790	4.XMFB-790	4.XMFS-790	4.0000-791
100 mg – 1 kg	25	2.110,9 g	4.XMFP-810	4.XMFB-810	4.XMFS-810	4.0000-811
100 mg – 2 kg	27	6.110,9 g	4.XMFP-820	4.XMFB-820	4.XMFS-820	4.0000-821
100 mg – 5 kg	28	11.110,9 g	4.XMFP-830	4.XMFB-830	4.XMFS-830	4.0000-831
100 mg – 10 kg	29	21.110,9 g	4.XMFP-840	4.XMFB-840	4.XMFS-840	4.0000-841

\*solange der Vorrat reicht